

# Tylosin Tartrate

(10E,12E)-(3R,4S,5S,6R,8R,14S,15R)-14-[[6-deoxy-2,3-di-O-methyl-B-D-allopyranosyl]oxymethyl]-5-[[3,6-dideoxy-4-O-(2,6-dideoxy-3-C-methyl- $\alpha$ -L-ribo-hexopyranosyl)-3-dimethylamino-B-D-glucopyranosyl]oxy]-6-formylmethyl-3-hydroxy-4,8,12-trimethyl-9-oxoheptadeca-10,12-dien-15-olide.  
Tylosin A (Tylosin)

916.10 CAS RN®: 1401-69-0; UNII: YEF4JXN031.

» Tylosin Tartrate is a tartrate of a mixture of macrolide antibiotic substances, or the mixture of such substances, produced by the growth of *Streptomyces fradiae*, or by any other means. Its potency is not less than 800  $\mu$ g of tylosin per mg, calculated on the dried basis.

**Packaging and storage**—Preserve in well-closed containers, protected from light, moisture, and excessive heat. Store at 25°, excursions permitted between 15° and 30°.

**Labeling**—Label it to indicate that it is for veterinary use only.

**USP REFERENCE STANDARDS (11)**—

[USP Tylosin RS](#)

[USP Tylosin Tartrate RS](#)

**Change to read:**

**Identification**—

**A:** ▲*Spectroscopic Identification Tests* (197), *Infrared Spectroscopy*: 197K▲ (CN 1-May-2020)

**B:** The retention time of the major peak for tylosin A in the chromatogram of the *Test solution* corresponds to that in the chromatogram of the *Standard solution*, as obtained in the test for *Content of tylosins*.

**C:** It meets the requirements of the test for *Tartrate* (191).

**PH (791):** between 5.0 and 7.2 in a solution prepared by dissolving 0.25 g in 10 mL of carbon dioxide-free water.

**LOSS ON DRYING (731)**—Dry about 1 g, accurately weighed, in vacuum at a pressure not exceeding 5 mm of mercury at 60° for 3 hours: it loses not more than 4.5% of its weight.

**RESIDUE ON IGNITION (281):** not more than 2.5%, the charred residue being moistened with 2 mL of nitric acid and 5 drops of sulfuric acid.

**Limit of tyramine**—In a 25-mL volumetric flask, dissolve 50.0 mg of tylosin in 5.0 mL of a 3.4 g per L solution of phosphoric acid. Add 1.0 mL of pyridine and 2.0 mL of a saturated solution of ninhydrin (about 40 g per L). Close the flask with aluminum foil, and heat in a water bath at 85° for 30 minutes. Cool the solution rapidly to room temperature, and dilute with water to volume. Mix, and measure immediately the absorbance (see [Ultraviolet-Visible Spectroscopy \(857\)](#)) of the solution at 570 nm against a blank solution prepared in a similar manner. The absorbance is not greater than that of a standard prepared at the same time and in the same manner using 5.0 mL of a 35 mg per L solution of tyramine in a 3.4 g per L solution of phosphoric acid. If intended for use in the manufacture of parenteral dosage forms, the absorbance is not greater than that of a standard prepared at the same time and in the same manner using 5.0 mL of a 15 mg per L solution of tyramine in a 3.4 g per L solution of phosphoric acid.

**Content of tylosins**—

**Mobile phase**—Prepare a mixture of filtered 200 g per L of sodium perchlorate, previously adjusted with 1 N hydrochloric acid to a pH of 2.5  $\pm$  0.1, and acetonitrile (60:40). Make adjustments if necessary (see *System Suitability* under [Chromatography \(621\)](#)).

**Standard solution**—Dissolve an accurately weighed quantity of [USP Tylosin RS](#) in a mixture of acetonitrile and water (1:1) to obtain a solution having a known concentration of about 0.2 mg per mL. [NOTE—Prepare the *Standard solution* immediately before use.]

**Test solution**—Dissolve an accurately weighed quantity of Tylosin in a mixture of acetonitrile and water (1:1) to obtain a solution having a known concentration of about 0.2 mg per mL. [NOTE—Prepare the *Test solution* immediately before use.]

**Chromatographic system** (see [CHROMATOGRAPHY \(621\)](#))—The liquid chromatograph is equipped with a 290-nm detector and a 4.6-mm  $\times$  20-cm column that contains 5- $\mu$ m packing L1. The flow rate is about 1.0 mL per minute and the column temperature is maintained at 35°.

Chromatograph the *Standard solution*, and record the peak responses as directed for *Procedure*: the order of elution is tylosin C, tylosin B, tylosin D, and tylosin A with relative retention times of about 0.5, 0.6, 0.8, and 1.0 minutes, respectively; the resolution of the peaks representing tylosin D and tylosin A is not less than 2.0; the tailing factors are not more than 1.5; and the relative standard deviation for replicate injections is not more than 2.0%.

**Procedure**—Separately inject equal volumes (about 20  $\mu$ L) of the *Standard solution* and the *Test solution* into the chromatograph, record the

chromatograms over a period of time equivalent to 1.5 times the elution time of the main tylosin A peak, and measure the peak areas for all the peaks. Calculate the percentages of tylosin A, tylosin B, tylosin C, and tylosin D in the Tylosin taken by the formula:

$$100(r_i/r_s)$$

in which  $r_i$  is the area of the tylosin A peak, the tylosin B peak, the tylosin C peak, or the tylosin D peak, as appropriate, in the chromatogram obtained from the *Test solution*; and  $r_s$  is the sum of the areas of all the peaks in the chromatogram obtained from the *Test solution*: the content of tylosin A is not less than 80%; and the sum of the contents of tylosin A, tylosin B, tylosin C, and tylosin D is not less than 95%.

**Assay**—Proceed as directed for Tylosin under [Antibiotics—Microbial Assays \(81\)](#). Prepare the *Test Dilution* as follows. Transfer an accurately weighed quantity of Tylosin Tartrate, equivalent to about 250 mg of tylosin, to a 500-mL volumetric flask, add 50 mL of methanol, and swirl to dissolve. Dilute with *Buffer B.3* to volume, and mix. Transfer 4.0 mL of this solution to a second 500-mL volumetric flask, dilute with a mixture of *Buffer B.3* and methanol (1:1), and mix. This solution contains about 4 µg of tylosin per mL.

**Auxiliary Information** - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
TYLOSIN TARTRATE	<a href="#">Julie Zhang</a> Associate Science & Standards Liaison	BIO42020 Biologics Monographs 4 - Antibiotics
REFERENCE STANDARD SUPPORT	RS Technical Services <a href="mailto:RSTECH@usp.org">RSTECH@usp.org</a>	BIO42020 Biologics Monographs 4 - Antibiotics

**Chromatographic Database Information:** [Chromatographic Database](#)

**Most Recently Appeared In:**

Pharmacopeial Forum: Volume No. PF 31(5)

**Current DocID:** GUID-961B7FA3-3D99-4581-95AF-8BC6D90E36BC\_2\_en-US

**DOI:** [https://doi.org/10.31003/USPNF\\_M87010\\_02\\_01](https://doi.org/10.31003/USPNF_M87010_02_01)

**DOI ref:** [up49y](#)